

# Laboratory Mathematics

---

Sorry In Advance

# Lab Math Concepts

---

- Conversions
- Molecular Weight and Concentration
  - Normality, Equivalence
- % solutions
- Dilutions
- Specific gravity

# Conversions

---

Dimensional Analysis

# Conversion Factors

---

- How do we use them?
  - Eg. 80 mm is how many inches? (conversion factor is  $1'' = 2.54 \text{ cm}$ )
- This process is repeatable and defined

# Dimensional Analysis

---

- By using fundamental qualities and their relationships (***conversion factors***) we can do most problems involving conversions
  - E.g. Calculate the molarity of 8.00 g of NaCl in 125 mL of water

1.094 *M*

# Molecular Weight and Concentration

---

Moles/Molarity and Equivalents/Normality

# Molecular Weights and Concentrations

---

## Molarity

- Moles per liter of solution
  - A mole of substance is  $1 \times$  gram molecular weight (gmw)
  - 23 g of sodium in 1 L is 1 M
  - 40 g of calcium in 1 L is 1 M

## Normality

- Equivalents per liter of solution
  - An equivalent is gmw / valence (measure of reactive potential)
  - 23 g of sodium in 1 L is 1 N
  - 40 g of calcium in 1 L is 2 N
  - $40 \text{ gmw} / 2+ \text{ charge} = \text{EqWt} = 20\text{g}$

# Equivalent weight

---

- What is the equivalent weight of  $\text{H}_3\text{PO}_4$ ? How about mEq wt?

- Molar weight?  $(3 \times 1) + 31 + (16 \times 4) = 98 \text{ g / mol}$

$$\frac{98 \text{ g H}_3\text{PO}_4}{1 \text{ mol H}_3\text{PO}_4} \times \frac{1 \text{ mol H}_3\text{PO}_4}{3 \text{ Eq H}_3\text{PO}_4} =$$


$$= 32.7 \text{ mg/ mEq}$$

# Normality

---

- It's just another step: 350 g NaCl in 1.5 L of water is what normality?
- Do you miss working out molarity? You can do that first if you'd prefer:

# % Solutions

---

Or an antiquated way to write concentrations

# % Solutions

---

- % (W/W)- The number of grams of substance in 100 g of solution
  - Eg: 10%  $\text{Li}_2\text{CO}_3$  solution is 10g  $\text{Li}_2\text{CO}_3$  in 100 g of mixture
- % (V/V)- The number of mL in a total volume of 100 mL
  - Eg: 10% HCl solution is 10 mL HCl in 100 mL of solution
- % (W/V)- the number of g of substance in a solution of 100 mL volume
  - 10% NaCl solution is 10g of NaCl in 100 mL of solution

# Dilutions

---

Creating solutions of lower concentration

# Dilutions

---

- Expressed relative to total
  - 1:2
  - $\frac{1}{2}$
- Others may express sample to diluent



Both of these  
are X 2 dilutions  
(D2 in chem)

# Dilutions

---

- Combining sample and diluent allows us to use the following equation:

$$V_1C_1 = V_2C_2$$

- Example: 100  $\mu\text{L}$  of serum is combined with 400  $\mu\text{L}$  of saline. What dilution have we made?
  - Part / Total  $\rightarrow$  100  $\mu\text{L}$  / 500  $\mu\text{L}$   $\rightarrow$  1/5
- What is our dilution factor?
  - x5

# Dilutions Application

---

- We can only measure so high, but some people are sicker than that, imagine... what nerve!?
- Analyzer can measure lipase 3-300 U/L
  - Returns result of: 5000 > Linearity – INVALID RESULT
- What dilution should we make?
  - 1/11 will give us ~450 U/L still too high
  - 1/101 will give us ~45 U/L, that works, but is large and introduces a lot of possible error
  - 1/21 will give us 238 U/L PERFECT

# Dilution Application

---

- Making a 1/21 aka D21 aka x21
- Choose a good total volume (protip: a multiple of 21)
- 210 $\mu$ L is good for the fixed volume pipettes we have
  - If only 1/21<sup>st</sup> of our dilution is patient sample that is 10 $\mu$ L
  - The rest (200 $\mu$ L) is diluent
  - 10 $\mu$ L patient + 200  $\mu$ L diluent = 210 $\mu$ L
- Dilution = 1/21
- Answer returned as 247 U/L
  - 247 U/L x 21 = 5187 U/L

# Specific Gravity

---

Density by another name

# Specific Gravity

---

- Measure of density
- A ratio of masses, represents 1 mL of given liquid solution (i.e. urine) in relation to 1 mL of H<sub>2</sub>O
- Urine typically has S.G. 1.005-1.030
  - Assuming we have urine with SG of 1.015
    - 1 mL of that urine has a mass of 1.015 g