




Quality Control and Statistics



Lab Statistics

- Medical Lab Scientists rely on data
 - How to interpret numbers
 - “Descriptive Statistics”
 - Mean Median and Mode
 - These numbers describe the overall data set
 - In Gaussian distribution these are all the same





Lab Statistics

- After we know what the data set is, we must know more about how it is spread out

- $\sigma = \sqrt{\frac{\sum(\bar{x} - x_n)^2}{n-1}}$

- $CV = \frac{\sigma}{\bar{x}} \times 100$ Measure of precision



Lab Statistics

- Not every data set will have “normal distribution” about the mean
 - It is important to see when data sets do not fit this curve




FIGURE 3-6 Examples of normal (Gaussian), skewed, and bimodal distributions. The type of statistical analysis that is performed to analyze the data depends on the distribution (shape).




Lab Statistics

- Much of what we know about how well methods work is based off of a “reference method”
 - The best known method produces a “known” or “true” value that is taken to be 100% accurate



Lab Statistics


- Correlation Coefficient (r)
 - Strength of relationship
 - Should be $>.98$
- Any deviation from the 1:1 line (reference method) is due to error
 - Random error- unpredictable but ever present
 - With increasing sample size the impact of random error is decreased



Lab Statistics


- Bias

- It can be difficult to see a bias if it is subtle or not present across measuring range
- This one is kind of obvious



Lab Statistics

- Systematic error- consistently moves measurements in one direction
 - Constant systematic error
 - Y-intercept constant +6
 - Proportional systemic error
 - Slope indicates underestimation
 - As value increases so does difference



Lab Statistics

- Why test the test?
 - Method evaluation and validation




FIGURE 3-10 A flowchart on the process of method selection, evaluation, and monitoring. (Adapted from Westgard JO, Quam E, Barry T. *Basic QC Practices: Training in Statistical Quality Control for Healthcare Laboratories*. Madison, WI: Westgard Quality Corp.; 1998.)



Lab Statistics

- How to choose a method
 - The ability of a method to diagnose disease (or exclude one) is what we care about when selecting a method
 - Sensitivity vs specificity
 - Trade-off, as sensitivity increases *generally* the specificity will decrease



Lab Statistics

- Understanding Sensitivity and specificity

	People who we <i>know</i> have a condition	People who we <i>know</i> do not have a condition
People who tested positive	True Positives	False positives
People who tested negative	False negatives	True Negatives



Lab Statistics

- Sensitivity: $\text{True positives} / (\text{true positives} + \text{false negatives})$
- Specificity: $\text{True negatives} / (\text{true negatives} + \text{false positives})$
- Next level of statistics PPV and NPV
 - PPV- if you test positive what is the likelihood you have a condition
 - NPV- if you do not test positive how sure are we that you don't have disease



Lab Statistics

- $PPV = \text{True positives} / (\text{True} + \text{False positives})$
- $NPV = \text{True negative} / (\text{True} + \text{False negatives})$



Lab Statistics

- Evaluate this method for spec. sens. PPV and NPV

	People having heart attacks	People not having heart attacks
cTnT > 0.030 ng/mL	254	15
cTnT < 0.030 ng/mL	7	812

Sensitivity: 97%

Specificity: 98%

PPV: 94%

NPV: 99%



Quality Control or QC

- Reliability (what do we want in a method)
 - Repeatable
 - Accurate
 - Precise
 - Trustworthy
 - Consistent



Quality

- Every lab must establish Quality Assurance Plan (QAP)
 - This must detect, control, and prevent occurrence of errors
- Quality Assurance is more than QC
 - Covers preanalytic, analytic, and post analytic phases



Quality

- Preanalytic controls
 - Physicians must order
 - Patient prep. i.e. fasting, 24 hr collection etc
 - Patient ID
 - Specimen acquisition, transport, processing and storage



Quality

- Analytic controls
 - What methodologies are chosen
 - Standardization and calibration techniques
 - Documentation of linearity and variability
 - Monitoring instrumentation, reagents
 - This is where QC comes in





Quality

- Postanalytic controls
 - Notify of critical values
 - HIPAA maintained
 - Maintain records of values





Quality Control

- Regular monitoring of tests
 - Analyze controls along with patient specimens
 - Internal QC
 - Monitor linearity
 - Ensure test limits met
 - External QC
 - Prof. Testing
 - QC samples




Quality Control

- Standard Deviation
 - Material assayed for 20 days
 - Establish mean and SD



Quality Control

- Accuracy and Precision
 - Accuracy: Closeness of agreement between measured value of analyte and its TRUE value
 - Precision: Ability to produce same value for replicate measurements



Quality Control

- Materials to measure accuracy:
 - Calibrator
 - Standards
 - Assayed Material
 - Reestablish the true values
- Materials to measure precision:
 - Unassayed QC
 - Patients from other instruments



05917824160V31.0

TSH CalCheck 5

ACCURACY
cobas®

REF 05917824 160

CalCheck 5 Lot: 183507

Reagent Pack Lot 120341 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 124433 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 128649 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 137811 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 143583 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 159557 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 188368 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Reagent Pack Lot 190074 for Elecsys 2010, MODULAR ANALYTICS E170, cobas e 411, cobas e 601 and cobas e 602 analyzers

Level	Value	Range	Unit
Check 1	≤ 0.010	≤ 0.010	μIU/mL
Check 2	1.97	1.62 - 2.32	μIU/mL
Check 3	48.0	40.8 - 55.2	μIU/mL
Check 4	81.0	68.9 - 93.2	μIU/mL
Check 5	> 100	90.1 - > 100	μIU/mL

COBAS, COBAS E and ELECSYS are trademarks of Roche.


© 2016, Roche Diagnostics

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.comDistribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

PRECISION

CHEM	TBIL	ALB	BUN	CA	CHOL	CL	CO2	CREAT	GLU	HDL	K	MG	NA	PHOS	TP	TRIG	TECH ID
LIMITS ±	0.2 or 20%	0.4 or 10%	2 or 9%	0.4 or 5%	4 or 10%	2 or 5%	3 or 15%	0.04 or 10%	6 or 10%	4 or 15%	0.5	0.2 or 10%	4	0.2 or 10%	0.2 or 5%	4 or 5%	
N756A/N1			20	7.6		103	24	1.50	103	31	7.1		144		6.2		cn
N756B/N2	0.5	3.1	22	7.7	114	102	23	1.52	107		7.1	2.0	144	3.4		105	cn
N766A/N3			21	7.7		104	23	1.49	106	30	7.2		146		6.0		cn
N766B/N4	0.5	3.2	21	7.4	110	104	23	1.49	107		7.2	2.0	145	3.4		104	cn
S756A/S1			20	7.7		104	22	1.45	105	30	7.2		145		6.0		cn
S756B/S2	0.5	3.1	21	7.6	110	105	22	1.49	105		7.1	1.9	146	3.4		98/101	cn
S766A/S3			20	7.6		104	22	1.52	106	30	7.1		145			5.8/5.9	cn
S766B/S4	0.5	3.2	21	7.6	112	104	22	1.53	105		7.2	1.9	145	3.3		100	cn
VAR.	0	0.1	2	0.3	4	≤ 5%	2	4.0%	4	1	0.1	0.1	2	0.1	6%	5.2%	3.9%

PRECISION



QC values

- Abnormal and normal values
- Range representative of assay
- Calculating the SDs
 - $\sigma = \sqrt{\frac{\sum(\bar{x} - x_n)^2}{n-1}}$
 - $CV = \frac{\sigma}{\bar{x}} \times 100$ (a universal SD language)
- Measures of precision



QC Values

- How is CV useful?

- Example:

- Potassium QC

Level 1: $\bar{x} = 2.6$ $\sigma = 0.1$ $CV = 3.85\%$

Level 3: $\bar{x} = 7.7$ $\sigma = 0.2$ $CV = 2.60\%$

- Cholesterol QC

Level 1: $\bar{x} = 105$ $\sigma = 5$ $CV = 4.76\%$

Level 3: $\bar{x} = 273$ $\sigma = 9$ $CV = 3.30\%$



Practice

- Let's evaluate the precision of 2 methods:

	Method 1	Method 2
\bar{x}	100 units	138 units
σ	25 units	31 units
CV?	25	23

Less precise

More precise



Control Situations

- What does it take to be out of control?
 - Any value greater than ± 3 SD
 - Both Levels greater than ± 2 SD
 - One Control $< \pm 2$ SD and another between ± 2 -3 SD
 - This is acceptable once every 20 days under LJ rules
 - Trends: In/decreasing for 6 consecutive days on both levels
 - Shift: All controls on one side of mean for 6 consecutive days on both levels



Shifts & Trends

- Why shift?
 - Improper storage of standard
 - New standard prepared incorrectly
 - Reagent shifted to new level of activity (enzymes at a different temperature)
 - Timer is not working properly
 - Different lot of reagent or standard



Shifts & Trends

- Why Trend?
 - Gradual deterioration of standard or reagent
 - Continued temp. instability (evap, damage)
 - Protein accumulation



QC

- What is it good for?
 - Alerts tech a problem exists
 - Does not ID problem
 - Does not fix problem
 - Troubleshooting *should* ID and fix the problem



Other QC

- Cumulative Summary Control Chart
 - Used in coagulation-hematology
- Moving Averages
 - Also used in heme, how stable the average for patient analytes is, if it is moving too much it's a problem
- Westgard Multi-Rule Control Chart
 - The Modern interpretation of LJ charts



Westgard Rules

- 1_{2S} -Warning, 1 control >2 SD
- 1_{3S} - Fail 1 control >3 SD
- 2_{2S} - Fail 2 consecutive controls > 2 SD
 - Also if *both levels* are >2 SD once
- R_{4S} - Range between levels greater than 4 SD
- 4_{1S} -Fail when 4 consecutive exceed same SD
- 10_x -Fail when 10 consecutive are on same side of mean

