



# Cytochemical Stains

Pinal Patel, MLS (ASCP)<sup>CM</sup>


# Leukocyte Alkaline Phosphate (LAP)

- Distinguish the cells of leukemoid reactions with increase activity from CML with decreased activity
- The LAP enzyme is in the tertiary granules of segmented neutrophils, bands and metamyelocytes
- The LAP within the neutrophil hydrolyzes the substrate (naphthol-AS-BI phosphate)
- hydrolyzed substrate then couples with the dye (fast blue BB salt) and precipitates at the site of activity
- ↑staining ↑enzyme activity







**Negative LAP reaction**



**Positive LAP reaction**


# Myeloperoxidase (MPO)

- Peroxidase enzyme is present in the primary granules of neutrophils(beginning promyelocyte stage and continuing throughout maturation),eosinophils and to certain extent monocytes
- Red cells and lymphocytes do not exhibit MPO activity
- Useful in differentiating myeloblasts from lymphoblasts



# Sudan Black B (SBB)




- Distinguish acute myelocytic and monocytic leukemia from acute lymphocytic leukemia
- SBB dye is fat soluble, stains fat particles present in the primary and secondary granules of myelocytic and monocytic cells.
- Lymphoid cells are unable to stain
- The reaction increases with the maturity of the granulocytic cell



# Esterases

- Differentiate myelocytic and monocytic leukemia
- WBCs contain esterases, a group of lysosomal enzymes
- Two substrate esters
  - $\alpha$ -naphthyl acetate and  $\alpha$ -naphthyl butyrate
    - both nonspecific and strong activity with cells of monocytic origin
  - Naphthol ASD chloroacetate
    - specific and positive reaction with granulocytic cells





# Tartrate Resistant Acid Phosphatase- TRAP

- Acid phosphatase is lysosomal enzyme present in many tissues and hematopoietic cells
- Seven isoenzymes have been classified as nonerythrocytic
- Isoenzyme 5 is tartrate resistant, produced in abundance by hairy cells



## Acute Leukemia Cytochemical Reaction Chart

---

Condition	MPO	SBB	NASDA	ANBE	ANAE
ALL	-	-	-	-/+ (focal)	-/+ (focal)
AML	+	+	+	-	-
AMML	+	+	+	+ (diffuse)	+ (diffuse)
AMoL	-	-/+	-	+ (diffuse)	+ (diffuse)
Megakaryocytic leukemia	-	-	-	-	+ (localized)