

# ABO Blood Group System

Andrea Nadas MBA, MLS  
(ASCP)<sup>CM</sup>




**Cleveland Clinic**













# Discovery of ABO Blood Group System

- 1901 by Karl Landsteiner
- Mixed cells + serum of associates
- Performed first forward and reverse grouping




# ABO Antigens and Antibodies


	Group A	Group B	Group AB	Group O
Red blood cell type				
Antibodies in Plasma	 Anti-B	 Anti-A	None	 Anti-A and Anti-B
Antigens in Red Blood Cell	 A antigen	 B antigen	 A and B antigens	None



# ABO Compatibility for RBCs



# ABO Compatibility for Plasma



# Forward Grouping (Front Type)



1 drop reagent  
(Anti-A, B, or  
D)



1 drop 3-5%  
patient red cell  
suspension













Centrifuge 30  
seconds




# Forwarding Grouping Results

Blood Type	Anti-A + Patient RBCs	Anti-B + Patient RBCs
O	0	0
A	4+	0
B	0	4+
AB	4+	4+

	Group A	Group B	Group AB	Group O
Red blood cell type				
Antibodies in Plasma	 Anti-B	 Anti-A	None	 Anti-A and Anti-B
Antigens in Red Blood Cell	 A antigen	 B antigen	 A and B antigens	None



# Reverse Grouping (Back Type)



1 drop A1 or B  
cells



2 drops  
patient plasma













Centrifuge 30  
seconds



# Reverse Grouping Results

Blood Type	A1 cells + Serum	B cells + Serum
O	4+	4+
A	0	4+
B	4+	0
AB	0	0

	Group A	Group B	Group AB	Group O
Red blood cell type				
Antibodies in Plasma	 Anti-B	 Anti-A	None	 Anti-A and Anti-B
Antigens in Red Blood Cell	 A antigen	 B antigen	 A and B antigens	None





# Full ABO Blood Type Results

Blood Type	Forward Grouping		Reverse Grouping	
	Anti-A	Anti-B	A1 cells	B cells
O	0	0	4+	4+
A	4+	0	0	4+
B	0	4+	4+	0
AB	4+	4+	0	0



# Frequency of ABO Blood Types in the Population

Most  
Common


Blood Group	Whites	African American	Hispanic	Asian
O	45%	50%	56%	40%
A	40%	26%	31%	28%
B	11%	20%	10%	25%
AB	4%	4%	3%	7%

Least  
Common




# ABO Antibody Characteristics

	ABO
Enzymes	Increased
IgM vs. IgG	Anti-A and Anti-B = IgM (Anti-A,B = IgG)
Cold or Warm Reacting	4°C (some 37)
Natural vs. Immune	Natural
Hemolytic Transfusion Reaction (HTR)	Yes
Hemolytic Disease of the Fetus and Newborn (HDFN)	Yes
Binds Complement	Yes
Dosage	No




# ABO Antibody Levels Throughout Life




# ABO Antigen Characteristics

- Develop in early fetal life
- Newborns = 20-50% of antigen sites
- A and B expression fully developed at 2-4 years
- Formed on:
  - RBC membranes
  - Endothelial cells
  - Platelets
  - Lymphocytes
  - Epithelial cells



# Inheritance of ABO Genes

- Genes: *A*, *B*, or *O*
- Chromosome 9
- *O* = amorph
- Codominant Expression



# ABO Genotype vs. Phenotype

Genotype	AA	AO	BB	BO	AB	OO
Phenotype	A		B		AB	O

Example: same phenotypes, different genotypes:

	A		
O		A	A
	O	AO	AO
	O	AO	AO


	A		
O		A	O
	O	AO	OO
	O	AO	OO




# Inheritance of Hh and Se genes

- FUT1 (H gene)
- FUT2 (Se gene)
- Closely linked on Chromosome 19
- Code for glycosyltransferases

TYPE 1 Secretions



TYPE 2 RBCs






# Formation of “O” Blood Type (RBCs)

Genotype: *HH* or *Hh*, *OO*

Gene	Glycosyltransferase	Immunodominant sugar	Antigen
<i>H</i> ( <i>FUT1</i> )	$\alpha$ -2-L-fucosyltransferase	L-fucose	H


Type O



# Formation of “A” Blood Type (RBCs)

Genotype: *HH* or *Hh*, *AA* or *AO*


Gene	Glycosyltransferase	Immunodominant sugar	Antigen
<i>H</i> ( <i>FUT1</i> )	$\alpha$ -2-L-fucosyltransferase	L-fucose	H
<i>A</i>	$\alpha$ -3-N-acetylgalactosaminyltransferase	N-acetyl-D-galactosamine (GalNAc)	A



# Formation of “B” Blood Type (RBCs)

Genotype: *HH* or *Hh*, *BB* or *BO*


Gene	Glycosyltransferase	Immunodominant sugar	Antigen
<i>H</i> ( <i>FUT1</i> )	$\alpha$ -2-L-fucosyltransferase	L-fucose	H
<i>B</i>	$\alpha$ -3-D-galactosyltransferase	D-galactose	B



# Formation of “AB” Blood Type (RBCs)

Genotype: *HH* or *Hh*, *AB*

Gene	Glycosyltransferase	Immunodominant sugar	Antigen
<i>H</i> ( <i>FUT1</i> )	$\alpha$ -2-L-fucosyltransferase	L-fucose	H
<i>A</i>	$\alpha$ -3-N-acetylgalactosaminyltransferase	N-acetyl-D-galactosamine (GalNAc)	A
<i>B</i>	$\alpha$ -3-D-galactosyltransferase	D-galactose	B



A antigen sites: 600,000

B antigen sites: 720,000




# Formation of Blood Type in Secretions

- Glycoproteins
- Inherit SeSe or Sese

80% Secretors (SeSe/Sese)  
20% Nonsecretors (sese)

Gene	Glycosyltransferase	Immunodominant sugar	Antigen
Se ( <i>FUT2</i> )	$\alpha$ -2-L-fucosyltransferase	L-fucose	Soluble H
A or B	$\alpha$ -3-N-acetylgalactosaminyltransferase $\alpha$ -3-D-galactosyltransferase	N-acetyl-D-galactosamine (GalNAc) D-galactose	Soluble A or B



# H Antigen Present by Blood Type

$O > A2 > B > A2B > A1 > A1B$

Greatest  
Amount of H

Least Amount  
of H



# Anti-H in A<sub>1</sub> and A<sub>1</sub>B Individuals

- Most H antigen is converted to A and B
- H is so well hidden patients can make anti-H
- Reacts best with O cells
- Interferes with antibody screens due to type O reagents

Anti-H (in A <sub>1</sub> and A <sub>1</sub> B)
IgM
Reacts at 4°C to room temp
Naturally occurring
No HTR - Insignificant
No HDFN



# A Subgroups



	A <sub>1</sub>	A <sub>2</sub>
Reactivity	Anti-A and Anti-A <sub>1</sub>	Anti-A only
Frequency	80%	20%
Inheritance	<i>A</i> gene normal	<i>A</i> gene mutated (single base substitution and single base deletion)
Enzyme	α-3-N-acetylgalactosaminyl-transferase	α-3-N-acetylgalactosaminyl-transferase with altered active site (less effective)
Immunodominant sugar	N-acetyl-D-galactosamine (GalNAc)	N-acetyl-D-galactosamine (GalNAc)
Produce Anti-A <sub>1</sub>	No	1-8% of A <sub>2</sub> individuals 22-35% of A <sub>2</sub> B individuals





# A<sub>1</sub> vs A<sub>2</sub> Structure

- 4 different H antigens:
  - H<sub>1</sub> and H<sub>2</sub> = unbranched straight chains
  - H<sub>3</sub> and H<sub>4</sub> = complex branched chains
- Terminal sugars are identical
- H<sub>1</sub> and H<sub>2</sub> converted to A<sup>a</sup> and A<sup>b</sup> antigens by both A<sub>1</sub> and A<sub>2</sub> enzymes (A<sub>2</sub> is less efficient)
- H<sub>3</sub> and H<sub>4</sub> converted to A<sup>c</sup> and A<sup>d</sup> antigens by A<sub>1</sub> enzyme and very poorly by A<sub>2</sub> enzymes or not at all



# $A_1$ and $A_2$ Quantitative and Qualitative Differences

---

## Quantitative

$A_2$  has decreased # of antigen sites ( $A_1 = 1,170,000/\text{RBC}$ ,  $A_2 = 290,000/\text{RBC}$ )

$A_2$  has decreased amount of transferase enzymes

$A_2$  has decreased amount of branched A antigens

## Qualitative

Differences in precursor chains

Differences in transferase

$A_2$  individuals can form an anti- $A_1$



# Anti-A<sub>1</sub> Lectin Testing






- Anti-A<sub>1</sub> lectin: reagent made from seed extract of *Dolichos biflorus* plant
  - Agglutinates A<sub>1</sub> but not A<sub>2</sub> cells
- Reactions of Patients' RBCs with:

	Anti-A (Anti-A plus Anti-A <sub>1</sub> )	Anti-A <sub>1</sub> lectin
A <sub>1</sub>	+	+
A <sub>2</sub>	+	0



# The Bombay Phenotype ( $O_h$ )


- Results from double dose of  $h$  gene ( $hh$ ) – extremely rare
  - Prevalence in India = 0.048%
  - Prevalence in Europe = 0.0003%
- No  $H$  gene means no H antigens which means no ABO genes expressed
- No H antigen also means they form anti-H

O	$O_h$
	
H only	None
	
Anti-A, Anti-B, Anti-A,B	 +Anti-H



# Bombay Phenotype Genetics

- Mutation in *FUT1* (*H* gene)
  - Silenced – unable to code for fucosyltransferase
  - **No**  $\alpha$ -2-L-fucosyltransferase, **no** L-fucose attached, **no** H substance present
- FUT2 gene (*Se* gene) also silenced
  - No ABO antigens form in secretions
- Normal ABO genes, but can't form ABO antigens due to lack of H
  - Genes written as  $Oh$ ,  $Oh^A$ ,  $Oh^B$ ,  $Oh^{AB}$



# Bombay Phenotype Testing

- Phenotype as O blood group (lack all ABH antigens)
- Screen is all positive
- Use Anti-H lectin (extract of *Ulex europaeus*)
  - Reacts with type O, but not with Bombay

Blood Group	Forward Group			Reverse Group		
	Anti-A	Anti-B	Anti-A,B	A1 Cells	B Cells	O Cells
A	+	-	+	-	+	-
B	-	+	+	+	-	-
AB	+	+	+	-	-	-
O	-	-	-	+	+	-
Bombay	-	-	-	+	+	+



# Anti-H in Bombay

	Anti-H (in A1 and A1B)	Anti-H in Bombay
IgM/IgG	IgM	IgM
Temperature of Reactivity	4°C and Room Temperature	Strongly at 37°C
Natural or Immune	Natural	Natural
Transfusion Reactions	No – Insignificant	Yes – binds complement & causes RBC lysis
HDFN	No	No

\*Bombay patients must receive Bombay blood (H negative)



# Para-Bombay

	Gene	RBC Antigens	Antigens in Secretions
Bombay	<b><i>hh sese</i></b> Silenced <i>FUT1</i> and <i>FUT2</i>	None	None
Parabombay (Red cell H partially deficient, nonsecretor)	<b><i>hh (weak variant) sese</i></b> Mutated <i>FUT1</i> gene With or without active <i>FUT2</i> gene	Very small amounts of H, A, and B	None
Parabombay (Red cell H deficient, secretor)	<b><i>hh Se</i></b> Silenced <i>FUT1</i> ( <i>H</i> ) gene Active <i>FUT2</i> ( <i>Se</i> ) gene	Little to none of H, A, and B (absorbed onto RBC from plasma)	H, A, and B







# THE FUTURE OF HEALTHCARE SINCE 1921

