





CLINICAL ENZYMES

Ryan Collison, MLS(ASCP)^{CM}SC^{CM}



Amylase (AMY/AMYL)

- Hydrolase that breaks down starch and glycogen
 - Requires Cl^- and Ca^{2+} activators
 - Acinar cells of pancreas and salivary glands are sources of amylase
 - Smallest enzyme, found in urine
 - Salivary amylase deactivated by stomach
 - Acute pancreatitis- AMY rises in 5-8 hours, peaks @ 24, and normal at 3-5 days



Amylase


- Somogyi Units
 - *Specifically amount of reducing sugar liberated (multiple products) We don't use.*
 - *U/L with our own ref. range*
- Measurement Methods
 - *Separate out salivary with wheat germ lectin*
 - *Immunoassay for specific Isoenzymes*
- Macroamylasemia: benign

| TABLE 13-6 AMYLASE METHODOLOGIES | |
|----------------------------------|---|
| Amylolytic | Measures the disappearance of starch substrate |
| Saccharogenic | Measures the appearance of the product |
| Chromogenic | Measures the increasing color from production of product coupled with a chromogenic dye |
| Continuous monitoring | Coupling of several enzyme systems to monitor amylase activity |

Copyright © 2014 Wolters Kluwer Health | Lippincott Williams & Wilkins


Lipase (LIPA/LPS)

- Hydrolyzes ester bonds to produce alcohol and fatty acids
 - *Specifically partially hydrolyzes triglyceride into 2-monoglyceride and 2 fatty acids*
 - *Pancreatic lipase is specific for the fatty acids at positions 1 & 2*
 - Substrate must be emulsified to occur, bile salts and colipase accelerate



Lipase

- Lipase is primarily found in pancreas
 - *Specific for pancreas!*
 - *Increases 4-8 hours after acute pancreatitis*
 - *Peaks at 24 hours*
 - *Normal at 8-14 days*
 - Longer lasting marker for pancreatitis
- Assay
 - *Early method Cherry-Crandall used olive oil to measure*
 - *Turbidimetric assay sees decrease in turbidity caused by fats*
 - *Colorimetric also based on glycerol kinase coupled reaction*




Pancreatitis

- Exocrine secretions digest foods
 - *Safeguards prevent activation of enzymes*
 - *Intrapancreatic secretion and activation*
 - “autodigestion”
- Precipitated by alcohol abuse, gallstones, trauma, [extremely] high triglycerides
 - *Treatments usually supportive*
 - *TPN only has given way to EN (Enteral Nutrition)*



Alkaline phosphatase (ALP)

- Group of enzymes that catalyze hydrolysis of phosphomonoesters at alkaline pH
 - *Liberates inorganic phosphate from organic molecule*
 - *Production of alcohol as a result*



- *Optimal pH 9-10 but varies with substrate*
 - Requires Mg^{2+} as activator
- Present on most cell outer surfaces
 - *Liver, bone, spleen, intestine, placenta and kidneys are highest*




Alkaline phosphatase

- Diagnostic Significance (specificity?)
 - *Liver and Bone*
 - Biliary tract obstruction- inducible enzyme will rise high
 - Bone disorders
 - *Paget's, osteomalacia & rickets, hyperparathyroidism, osteogenic sarcoma*
 - *Also present after fracture and during bone growth*
 - *Pregnancy- rises significantly from placental origin*



Alkaline phosphatase


- *Isoenzymes*
 - Electrophoresis to separate, Liver fastest, then bone, placental, intestinal
 - Heat stability, measure before and after 56°C for 10'
 - *Liver lives if <20% left, bone is source if >20% liver is source*
 - Placental even more, will resist 65°C for 30'
 - Chemical inhibition
 - *Phenylalanine inhibits intestinal and placental ALP more than liver and bone*
 - Regan and Nago isoenzymes- carcinoplacental ALP
 - *Occur in 3-15% of cancer patients, ectopically produced*
- Assay: leverage nonspecificity
 - *Bower & McComb*



Acid Phosphatase (ACP)


- Hydrolase that catalyzes similar reactions to ALP

- *Ideal pH is approx. 5.0*



- *Found in prostate, bone, liver, spleen, kidney, RBCs, platelets*
 - Prostate is largest source
 - Insensitive- Only detects prostate cancer after metastasis
 - *Replaced by PSA*
 - Inhibition: Prostatic ACP inhibited by tartrate
 - *Other usage? Rape kit test for ACP + up to 4 days*

Phosphatases Side-by-Side



Alkaline Phosphatase

- Optimal pH 9-10
- Tissue source: Liver, bone, spleen, intestine, placenta, kidneys
- Undemanding preanalytics
- Routine part of CMP, HFP

Acid Phosphatase


- Optimal pH ~5
- Tissue source: prostate
- Special storage (buffered, frozen, NO HEMOLYSIS)
- Esoteric test, rarely performed

Aspartate aminotransferase (AST)

- Transferase/transaminase


- Moves amino group from aspartate to α -keto acid
 - May be called SGOT (old name)
- Requires pyridoxal phosphate coenzyme
 - Ketoacids formed are eventually used in TCA cycle
- Mostly found in heart, liver, skeletal muscle, with little in kidney, pancreas, RBCs
 - Largely raised in liver disorders and skeletal muscle issues
 - Also seen in cases of pulmonary emboli
 - In cirrhosis some elevation 4xULN vs. hepatitis 100xULN

- Assay:




Alanine aminotransferase (alt)

- Another transferase similar to AST
 - *Also requires pyridoxal phosphate coenzyme*
 - *May also be called SGPT*
- Distributed widely by highly concentrated in liver
 - *More specific for hepatic disorders*
 - *ALT often rises higher than AST, remains high longer*




- Assay:

occurs best at pH 7.3-7.8



Transferases Side-by-Side



AST

- Found in liver, heart, skeletal muscle, RBCs
- Forms oxaloacetate and glutamate
- Requires pyridoxal phosphate coenzyme
- Rises in cases of pulmonary emboli, cirrhosis (mild), hepatitis (steep)

ALT

- Very liver specific
- Forms pyruvate and glutamate
- Requires pyridoxal phosphate coenzyme
- Rises higher than AST, but in only from liver disorders

γ -Glutamyltransferase (GGT)

- Transfers γ -glutamyl group from peptides to amino acids
 - *Clinically only important for the liver and biliary tract*
- Very sensitive to hepatobiliary disorders
 - *High elevations largely seen in blockage of biliary tract*
- Inducible enzyme
 - *Also raised by some drugs (warfarin, phenobarbital, phenytoin)*
 - *May indicate chronic alcoholism, monitor compliance*
- Assay conversion to a colored compound
 - *γ -glutamyl p-nitroaniline (GGPNA)*



Creatine Kinase (CK)

- ATP regeneration in muscle
- Dimeric enzyme composed of M and/or B subunits
 - *CK1 (BB)- Brain, most anodal*
 - Short half-life doesn't contribute to total
 - *CK2 (MB)- Heart*
 - Predominantly in heart, but also skeletal muscle
 - *CK3 (MM)- Muscle, most cathodal*
 - Predominantly skeletal, but also cardiac muscle






FIGURE 13-5 Electrophoretic migration pattern of normal and atypical creatine kinase (CK) isoenzymes.

Creatine Kinase


■ Significance

- *Frequently seen in disorders of cardiac and skeletal muscle*
 - MI, rhabdomyolysis, muscular dystrophy
 - *Extremely high in Duchenne MD 50-100x ULN*
- *Not very specific for MI, also seen in strokes, seizures, nerve degeneration, CNS shock etc...*
- *Separation into isoenzymes improves specificity*
 - In healthy person <6% of total is MB, mostly MM ~ 94%
 - CK-MB is high in cardiac tissue
 - *Begin to rise after 4-8 hours, peak at 12-24 hours, normal at 48-72 hours*
- Replaced by *troponin*, released earlier and lasts longer




Creatine Kinase

- Assay: can go either way with PK & LD




- Or with Hexokinase and G6-PDH



Lactate Dehydrogenase (LD)

- Converts lactic acid to pyruvic acid and generates NADH



- High amounts found in the heart, liver, muscle, kidney, and RBCs
 - *Elevated in numerous disorders*
 - Cardiac, hepatic, skeletal muscle, renal, hematological, and neoplasms
 - Highest total levels from anemias

Lactate dehydrogenase

- Viral hepatitis, cirrhosis
 - *Minimal increase*
- AMI & pulmonary infarct
 - *Same minimal increase*
 - LD rises within 12-24 hours, peaks at 24-48 hours, stays up for 10 days
- Leukemia
 - *Marked elevation in ALL possible*
- Isoenzymes allow greater specificity


| TABLE 13-5 | | LACTATE DEHYDROGENASE (LD) ISOENZYMES AS A PERCENTAGE OF TOTAL LD |
|---|--|---|
| ISOENZYME | | % |
| LD-1 | | 14–26 |
| LD-2 | | 29–39 |
| LD-3 | | 20–26 |
| LD-4 | | 8–16 |
| LD-5 | | 6–16 |
| Source: Lott JA, Stang JM. Serum enzymes and isoenzymes in the diagnosis and differential diagnosis of myocardial ischemia and necrosis. <i>Clin Chem</i> . 1980;26:1241. | | |

Lactate Dehydrogenase

- Each enzyme consists of 4 peptide chains
 - *LD-1 (HHHH)- Heart & RBCs*
 - *LD-2 (HHHM)- Heart and RBCs*
 - *LD-3 (HHMM)- Lung, WBCs, spleen, pancreas*
 - *LD-4 (HMMM)- Liver*
 - *LD-5 (MMMM)- Skeletal muscle*
- In AMI and intravascular hemolysis $LD-1 > LD-2$
 - *This is not the normal state, so this is known as the LD flipped pattern*
- LD-6?
 - *Alcohol dehydrogenase that migrates to LD-5 spot*

Glucose-6-Phosphate dehydrogenase (G-6-PD)

- Oxidoreductase that is part of the pentose phosphate shunt glucose metabolism
- Protects cells from oxidative damage
 - *Tissue sources: adrenal cortex, spleen, thymus, lymph nodes, mammary gland, and RBCs*
 - Little activity in normal serum
 - *Most of its importance is in the RBC*
 - Maintains NADPH in reduced form
 - *Insufficiency results in low NADPH, exposure to oxidizing agents, cells will burst from damage*
 - *Deficiency is Sex-Linked trait*
 - Drug induced hemolytic anemia



✱ Sample is hemolysate

Cytochrome Oxidase

- Cytochrome p450 oxidase (CYP 450)
 - *Superfamily of enzymes (absorb light at 450 nm)*
 - Metabolize more than 50% of all drugs
 - More than 500 different enzymes
 - *Different alleles allow different speeds of drug metabolism*
 - *Also found in steroid forming tissues*
 - Synthesize steroid hormones from cholesterol
- Pharmacogenomics depends upon these enzymes to determine effective therapies



Cholinesterases

- Break down choline neurotransmitters
- Target of organophosphate pesticides and some therapeutic drugs
 - *Acetylcholinesterase (AChE) intracellular*
 - *Pseudocholinesterase or Serumcholinesterase (SChE) extracellular*
 - Butyrylcholinesterase (BChE) preferred, new name
 - *Butyrylcholine a synthetic choline used to distinguish from AChE*
- Sarin, Tabun, Soman, VX
 - *Chemical weapons, nerve gas*
- Chronic toxicity- wide neurological symptoms
- Ref: 4-12k U/L
 - *Symptoms with a decrease of 40%*

