



Diagnostic Parasitology

Medically Important Amebae

Kaitlin Landfried, MLS(ASCP)^{CM}

Disclaimer

- This presentation was meant to provide students with both didactic and laboratory skills as they apply to clinical parasitology. It is meant for educational purposes only and does not represent Cleveland Clinic views or practices.
- The presentation contains images and other references copyrighted by another entity or person and credits shall be given to the rightful owners of the materials and I claim no copyright to the said content.
- Most of the information was adopted from the Textbook of Diagnostic Microbiology by Mahon & Lehman (see citation) but condensed for bite sized learning.

Protozoa

- Unicellular parasites
- Classified according to their motility organelles
 - Amebae move by pseudopodia
 - Ciliates move by cilia
 - Flagellates move by flagella
 - Sporozoa are nonmotile
- All reproduce asexually except for the sporozoa

Cyst: nonmotile form, resistance to environmental factors
Trophozoite: feeding, motile form, replicates in the host and responsible for causing damage

The Amebae

Intestinal

- Six medically important members
 - *Entamoeba histolytica/dispar/moshkovkii*
 - *Entamoeba hartmanni*
 - *Entamoeba coli*
 - *Endolimax nana*
 - *Iodamoeba butschlii*
 - *Blastocystis* species
- Includes pathogens and nonpathogens
- Humans are the only known host except for *Blastocystis* species
- Cyst is the infective form
- Fecal oral transmission

Tissue


- Three medically important members
 - *Naegleria fowleri*
 - *Acanthamoeba* species
 - *Balamuthia mandrillaris*
- Free-living in soil and water
- Trophozoite is the infective form

Species identification is determined based on size, number of nuclei, nuclear structure, and presence of specific internal structures

Amebae Microscopic Identification


STEP ONE

- Cyst vs Trophozoite
 - Amebic cysts can have more than one nucleus, depending on maturity, and lack motility organelles
 - Amebic trophozoites have one nucleus and motility organelles




Amebae Microscopic Identification


- Cyst or trophozoite?




Cyst



Trophozoite



Trophozoite




Cyst

Amebae Microscopic Identification


STEP TWO

- Describe the **karyosome**: a chromatin mass within the nucleus of protozoa
 - Shape: regular(circular) vs irregular
 - Size: dot-like (small) or blot-like (large)
 - Location: central or eccentric




Amebae Microscopic Identification

- Describe the karyosome.




regular
dot-like
centrally located



irregular
blot-like
centrally located



regular
blot-like
centrally located




regular
dot-like
eccentrically located

Amebae Microscopic Identification


STEP THREE

- Describe the **peripheral chromatin**: **DNA present on the nuclear membrane of some protozoa.**
 - Present/absent
 - Distribution: even vs uneven




Amebae Microscopic Identification


- Describe the peripheral chromatin




present
evenly distributed



absent



absent




present
unevenly distributed

Amebae Microscopic Identification


STEP FOUR

- Describe unique features:
 - Chromatoidal bars: condensed RNA found in the cytoplasm of some amebic cysts
 - a bar-like structure rounded or splintered
 - Vacuoles
 - Ingested red blood cells




Amebae Microscopic Identification


- Describe the unique features




ingested RBCs



rounded chromatoidal bar



splintered chromatoidal bar



large vacuole

Amebae Microscopic Identification

STEP FIVE


- Measure the size
 - Use the ocular micrometer to measure the size of the organism in μm
 - Some ameba look so similar the only way to differentiate is size
 - Always double check the size with the expected size for the parasite of suspicion

Amebae Microscopic Identification


Calibration of the Ocular Micrometer

Instructions

1. Insert the ocular micrometer into a 10X eyepiece. The ocular micrometer is divided into ocular divisions (OD).
2. Place the calibrated stage micrometer slide on the stage and focus on the scale. The stage micrometer has a calibrated scale which is divided into 0.1 millimeters (mm) and 0.01 mm units.
3. Adjust the field so the 0 line of the ocular micrometer (OM) scale is exactly superimposed upon the 0.0 line of the stage micrometer (SM) scale.
4. Without moving the stage micrometer, locate the point as far to the extreme right as possible where any two lines are exactly superimposed upon each other.
5. Count the number of divisions (mm) on the stage micrometer between the 0.0 line and the superimposed line to the far right.
6. Count the number of ocular divisions on the ocular micrometer between the 0 line and the superimposed line to the far right.
7. Divide the distance determined in step 5 by the number of ocular divisions in step 6 and multiply by 1000 to give the ocular micrometer units in μm .
8. Repeat steps 3 through 7 for each objective on the microscope.
9. If at any time the ocular micrometer is moved to a different microscope or a new objective is added to the microscope, the calibration procedure must be completed again.



ocular micrometer




stage micrometer

7.
$$\frac{\text{Length seen on stage micrometer}}{\text{\# of divisions counted on ocular micrometer}} \times 1000 = \text{ } \mu\text{m}$$

Use me!

Amebae Microscopic Identification

- Calibrate the ocular micrometer on 100x



7. Length seen
on stage micrometer


of divisions counted
on ocular micrometer

$\times 1000 = \text{ } \mu\text{m}$

$$\frac{0.2\text{mm}}{27} \times 1,000 = 7.4 \mu\text{m}$$

Amebae Microscopic Identification

- Determine the size of the object in micrometers (viewed on 100x)



Ocular distance x calibration =
size of object in micrometers

$$15 \times 7.4\mu\text{m} = 111\mu\text{m}$$

Amebae Microscopic Identification

STEP ONE

Form (Eluffy)

STEP TWO

Karyosome (Kittens)

STEP THREE

Peripheral chromatin (Purr)

STEP FOUR

Unique characteristic (Uncontrollably)

STEP FIVE

^{EVERY}
Size (Sunday)




Intestinal Amebae: *Entamoeba histolytica*/*dispar*/*moshkovkii*

- Three species are morphologically identical organisms
 - *Entamoeba histolytica* is considered a pathogen
 - *Entamoeba dispar* and *Entamoeba moshkovkii* are considered nonpathogenic
- *Entamoeba histolytica* is found worldwide
 - Asymptomatic colonization
 - Amebic dysentery
 - Acute: vague abdominal symptoms
 - Chronic: formation of ameboma
 - Extraintestinal amebiasis
 - Organism erodes the intestines and enters the circulation

Intestinal Amebae: *Entamoeba histolytica*/*dispar*/*moshkovkii*

Nucleus: **regular dot-like centrally located karyosome** with **evenly distributed peripheral chromatin**



Trophozoite


Size: 15-25 μm

Unique: presence of ingested RBCs

*ingested RBCs definitively identify the organism as *Entamoeba histolytica* *



TARGET



Cyst

Size: 10-20 μm


Unique: may contain chromatoidal bars with rounded ends

Intestinal Amebae: *Entamoeba histolytica*

- Other tests for identification
 - Tests that detect *E. histolytica* antigen in stool provide evidence of current infection
 - More sensitive than a microscopic examination of feces
 - These tests use the EIA method with monoclonal antibodies to proteins (e.g. serine-rich antigen or the galactose/ N-acetylgalactose adhesion lectins)
 - Multiplex molecular assays for GI syndromic testing are also available, and *E. histolytica* is often included


Intestinal Amebae: *Entamoeba hartmanni*

Nucleus: **regular dot-like centrally located karyosome** with **evenly distributed peripheral chromatin**



Trophozoite

Size: $\leq 12 \mu\text{m}$



Cyst

Size: $\leq 10 \mu\text{m}$

Unique: may contain chromatoidal bars with rounded ends

Intestinal Amebae: *Entamoeba coli*

Nucleus: **regular dot-like eccentrically located karyosome** with **uneven peripheral chromatin**




Trophozoite
Size: 15-50 μ m




Cyst
Size: 15-25 μ m
Unique: may contain chromatoidal bars with splintered ends.
Up to eight nuclei.

Intestinal Amebae: *Endolimax nana*

Nucleus: **irregular blot-like centrally located karyosome** without peripheral chromatin



Trophozoite
Size: 5-12 μ m




Cyst
Size: 5-12 μ m




Intestinal Amebae: *Iodamoeba butschlii*

Nucleus: **regular blot-like centrally** located **karyosome** surrounded by achromatic granules without **peripheral chromatin**



Trophozoite
Size: 6-20µm




Cyst
Size: 6-15µm
Unique: large vacuole.
one nucleus.



Intestinal Amebae: *Blastocystis* species

- Not exactly an ameba, but instead closely related to brown algae
- Controversy over pathogenicity
 - Considered a pathogen when patient is symptomatic, and only parasite found



Vacuolar form

Size: 5-15 μ m


Unique: large central vacuole up to 6 nuclei pushed up against the plasma membrane



Amebae Ova and Parasite Examination (OVAP)



Iodine wet preparation



F:

K:

P:


U:

S: 15µm


Trichrome



Amebae Ova and Parasite Examination (OVAP)




Iodine wet preparation




Trichrome



Amebae Ova and Parasite Examination (OVAP)



Iodine wet preparation



F:

K:

P:


U:

S: 25µm


Trichrome



Amebae Ova and Parasite Examination (OVAP)



Iodine wet preparation




F:
K:
P:
U: N/A
S: 6µm


Trichrome



Amebae Ova and Parasite Examination (OVAP)



Iodine wet preparation



F:

K:

P:


U:

S: 12 μ m

Trichrome



Amebae Ova and Parasite Examination (OVAP)



Iodine wet preparation

Trichrome




Tissue Amebae: *Naegleria fowleri*


- Primary amebic meningoencephalitis (PAM) in immunocompetent & immunocompromised
 - Trophozoite gains access to CNS by inhalation into the upper respiratory tract followed by the penetration of the nasal mucosa, which allows them to travel along the olfactory nerve to the brain
 - Headache, fever, stiff neck, nausea
 - Rapidly fatal
 - CSF will have increased neutrophils and proteins with decreased glucose

Tissue Amebae: *Naegleria fowleri*

Nucleus: regular blot-like centrally located karyosome without peripheral chromatin



Cyst



Trophozoite
Size: 10-35µm




Transient flagellate

Tissue Amebae: *Naegleria fowleri*


- Tests for identification
 - a) Motile trophozoites can be seen in a wet mount of CSF
 - b) Trophozoites can be seen on a Wright stain of CSF
 - c) Can be cultured on nonnutrient agar overlayed with *Escherichia coli*




a



b



c




<https://www.youtube.com/watch?v=2r8iJXwJQas>

Tissue Ameba: *Acanthamoeba* species

- Granulomatous amoebic encephalitis (GAE) in immunocompromised hosts
 - Months to years
- Amebic keratitis in immunocompetent due to improper storage and disinfection procedures (ex. wearing contacts while swimming)
- Cutaneous infections in AIDs patients
 - Chronic nonhealing lesions

Tissue Ameba: *Acanthamoeba* species


Nucleus: regular blot-like centrally located karyosome without peripheral chromatin



Trophozoite

Size: 15-45 μ m

Unique: Single central nucleus with blunt pseudopods with spinelike projections




Cyst

Size: 15-20 μ m


Unique: Single central nucleus with doubled walls and a wrinkle appearance

Tissue Ameba: *Acanthamoeba* species

- Tests for identification
 - Motile trophozoites can be seen in a wet mount of clinical specimens
 - Trophozoites and cysts can be seen on a permanent smears
 - Can be cultured on nonnutrient agar overlayed with *Escherichia coli*



<https://www.youtube.com/watch?v=qBldfJEJMK0>



Tissue Ameba: *Balamuthia mandrillaris*

- Primarily found in soil
- GAE and cutaneous infections in immunocompetent and immunocompromised hosts.


Nucleus: regular blot-like centrally located karyosome without peripheral chromatin



Trophozoite
Size: 30-60µm
Unique: single central nucleus
with broad pseudopods



Cyst
Size: 10-30µm
Unique: single central nucleus
multiwalled



Citations

- Mahon, C. R., & Lehman, D. C. (2023). *Textbook of Diagnostic Microbiology* (7th ed., pp. 639-707). Elsevier.
- Centers for Disease Control and Prevention (2019, November 20). DPDx-Laboratory Identification of Parasites of Public Health Concern. Retrieved November 13, 2023, from <https://www.cdc.gov/dpdx/az.html>