

# Anaerobic Bacteriology

Susan M. Harrington, Ph.D.  
2022

# Anaerobic bacteria

- Dominate indigenous flora
- Commonly found in infection
- Special procedures required for collection and transport
- Most commonly overlooked bacterial infection
- Endogenous or exogenous sources

# Major anaerobes encountered clinically

(see flow chart in text)

- **Spore-forming GPR**
  - *Clostridium* spp.
- **Non-spore-forming GPR**
  - *Actinomyces* spp.
  - *Cutibacterium* (*Propionibacterium*)spp.
  - *Mobiluncus* spp.
  - *Lactobacillus* spp.
  - *Eubacterium* spp.
  - *Bifidobacterium* spp.
- **Gram-positive cocci**
  - *Peptostreptococcus* spp.
- **Gram-negative bacilli**
  - *Bacteroides fragilis* group
  - Other *Bacteroides* spp.
  - *Porphyromonas* spp.
  - *Prevotella* spp.
  - *Fusobacterium* spp.
  - *Bilophila* spp.
  - *Sutterella* spp.
- **Gram-negative cocci**
  - *Veillonella* spp.

# Habitat of anaerobes

- Widespread: soil, marshes, lake and river sediments, ocean, sewage, foods, animals
- Humans
  - oral cavity around teeth
  - URT
  - GI tract
  - Orifices of GU tract
  - skin

# Incidence of anaerobes as indigenous flora in humans

<b>Gram Negative Anaerobes</b>	<b>Skin</b>	<b>URT</b>	<b>Mouth</b>	<b>Intestine</b>	<b>GU tract</b>
<i>B. fragilis group</i>	0	0	0	2	+/-
<i>Bilophila</i>	0	0	+/-	1	+/-
<i>Fusobacterium</i>	0	1	2	1	1
<i>Porphyromonas</i>	0	+/-	1	1	+/-
<i>Prevotella</i>	0	1	2	1	1
<i>Sutterella</i>	0	0	0	1	0
Other GNR, GNC	0	1	1	1	1

# Incidence of anaerobes as indigenous flora in humans

<b>Gram Positive Anaerobes</b>	<b>Skin</b>	<b>URT</b>	<b>Mouth</b>	<b>Intestine</b>	<b>GU tract</b>
<i>Clostridium</i>	0	0	+/-	2	+/-
<i>Actinomyces</i>	0	1	2	+/-	+/-
<i>Bifidobacterium</i>	0	0	1	2	+/-
<i>Eubacterium</i>	+/-	+/-	1	2	+/-
<i>Lactobacillus</i>	0	0	1	1-2	1-2
<i>Cutibacterium</i>	2	1	+/-	+/-	+/-
<i>Peptostreptococcus</i>	1	2	2	2	2

# Events that lead to infection


## Endogenous:

- Dental procedures
- GI surgery/traumatic abdominal injury
- Genital tract surgery, IUD use
- Antibiotic use

## Exogenous:

- Human or animal bites
- Traumatic surface wounds
- Improperly prepared/stored food

# Anaerobic Infections (often polymicrobial)



*C. acnes* – PJI  
*C. difficile* - AAC

*Clostridium* spp, *B. fragilis* group,  
*Peptostreptococci*,  
*Fusobacterium* spp.

*Peptostreptococci*,  
*Bacteroides* spp.,  
*Clostridium* spp,  
*Prevotella* spp.,  
*Actinomyces israelii*  
(IUD)

*Clostridium* spp.

Koneman 1997; Fig 14-1

& Mahon Text of Diagnostic  
Micro 6<sup>th</sup> Ed

ANA GNB &  
*Clostridium* spp.

*Peptostreptococci*,  
*Fusobacterium* spp.,  
*Porphyromonas* spp.

*Porphyromonas*,  
*Fusobacterium*,  
*Actinomyces*,  
*Finegoldia magna*,  
*Peptostreptococci*, *B. fragilis* group


*B. fragilis* group,  
*Fusobacterium* spp.,  
*Clostridium* spp.,  
*Peptostreptococci*



# Specimen –

## Indications for Anaerobic Infection

- Infection close to mucosal surface (pyogenic)
- Foul odor
- Large quantity of gas
- Sulfur granules
- Black discoloration of tissue
  - Pigments of *Porphyromonas* spp. & *Prevotella* spp. may cause blackening



## Appropriate for Anaerobes (Deep Sites)

- Aspirates (by needle)
- Blood and bone marrow
- Bronchoscopy by protected specimen brush
- Transtracheal aspirate
- Nasal sinus aspirate
- Deep tissue (organs, placenta, surgically obtained)
- IUD for Actinomyces
- Suprapubic urine

## Inappropriate for Anaerobes

- Cervical, vaginal & urethral secretions
- Prostatic or seminal fluid
- BAL & bronchial washes
- NP swabs & sinus washings
- Endotracheal aspirates
- Sputum
- Surface wounds or tissues
- Skin scrapings
- Stool or rectal swabs
- Urine (voided or catheter)

# Relationship of bacteria to oxygen

- Obligate aerobes
  - *Micrococcus*, *Pseudomonas* spp.
  - Require O<sub>2</sub> as terminal electron acceptor; no fermentation
- Facultative anaerobes
  - *E. coli*, *S. aureus*
  - Grow under aerobic or anaerobic conditions
- Microaerophiles
  - *Campylobacter* – optimal growth in 5% O<sub>2</sub>; 10% CO<sub>2</sub>
  - Require O<sub>2</sub> as terminal electron acceptor, yet do not grow in aerobic incubator (21% O<sub>2</sub>); grow minimally under anaerobic conditions

# Anaerobes

- Strict obligate - cannot tolerate  $>0.5\%$   $O_2$ 
  - *Clostridium haemolyticum*, *C. novyi* type B
- Moderate obligate - can tolerate 2-8%  $O_2$ 
  - *Bacteroides fragilis*, *Prevotella*-*Porphyromonas* groups, *Fusobacterium nucleatum*, *C. perfringens*
- Aerotolerant anaerobe – show scant growth on agar in room air or 5-10%  $CO_2$ ; good growth under anaerobic conditions
  - *Clostridium* *carnis*, *C. histolyticum*, *C. tertium*
  - *C. acnes*, *Lactobacillus*

# Culture for Anaerobes

- To preserve strict anaerobes culture within 20 min.
- Swabs are poor specimens.
- PRAS transport media for swabs, fluids & tissues is preferred. E swab.
- If specimen is small put the whole thing in transport & process for aerobes & anaerobes.
- Refrigerated/frozen unacceptable




# Primary Anaerobic Media

- Media should be RT; pre-reduced.
- Anaerobic Blood Agar (CDC) or Brucella Blood Agar
  - Supports growth of almost all obligate and facultative anaerobes
  - Vit K, hemin, yeast extract
- Kanamycin-Vancomycin Laked Blood Agar (KVLB)
  - Selective isolation of *Bacteroides* and *Prevotella*
  - Pigmentation of *Prevotella*





# Primary Anaerobic Media

- Bacteroides Bile Esculin Agar (BBE)
  - Selective, differential for *B. fragilis* group
    - Bilophila wadsworthia grows
  - 20% bile salts & gentamicin
  - Bile resistant organisms grow
  - Hydrolysis of esculin – black/brown colony formation/halo
    - *B. vulgatus* is esculin neg



# Anaerobic Media

- Thioglycolate Broth
  - Provides enrichment
  - Non-selective
  - Thioglycolate & cystine – reducing agents
  - Low conc of agar
  - Resazurin indicator (pink if oxidized)






# Anaerobic Media

- Phenylethyl Alcohol Blood Agar (PEA)
  - Usually w/sheep blood
  - Inhibits facultative GNBs
  - Inhibits swarming
- Cycloserine cefoxitin fructose agar (CCFA)
  - 500 µg/ml cycloserine; 16 µg/ml cefoxitin
  - *C. difficile*; yellow colony
  - Horse manure odor
  - UV light fluorescence

CCFA



# Differential Media

- Egg Yolk Agar
  - Egg emulsion
  - Identification
    - lecithinase (opaque) (Nagler reaction)
    - lipase (sheen)
    - protease (clearing)



EYA - Lipase



EYA - Lecithinase


# Anaerobic incubation

- 5% H<sub>2</sub>, 5-10% CO<sub>2</sub>, 85-90% N<sub>2</sub>
- Techniques
  - Systems that generate gas from envelope (pack):
    - jar
    - bag
    - box
  - Anoxomat
    - Gas generated by tank connected to jar
- Methylene blue indicator strip:  
white = ANA



# Anaerobic Incubation

- Anaerobic chamber (glove box)



gloveless

# Incubation Times: Anaerobic Cultures

- Incubated for 5 days
- Media usually includes: CDC & BBE or LKV in ANO<sub>2</sub>
  - BAP O<sub>2</sub>; Thio for selected sites.
- Incubated for 48 hours before first examination.
- Extended incubation:
  - *C. acnes* (7 – 14 days)
  - Actinomyces (14 days)

# Identification of Anaerobes

## Preliminary Tests


- Gram stain
- Colony morphology/differential growth
- Pigment
- Fluorescence
- Catalase (15% H<sub>2</sub>O<sub>2</sub>)
- Spot indole
- Nitrate
- Aerotolerance



PURE COLONY ISOLATE

Aerotolerance test:  
Compare growth of  
aerobic and anaerobic  
plates:

BAP CO<sub>2</sub>  
CDC ANO<sub>2</sub>



<sup>x</sup>Optional but highly recommended  
RLBA = Rabbit laked blood agar

# Identification of anaerobes

- Special potency antibiotic disks
- Commercial biochemical systems
  - PRAS biochemicals (std method, exp, labor intensive)
  - API 20A (24-48 h)
    - Saccharolytic, rapid growing (*B. frag*, *Clostridium* spp.)
  - RapID ANA II (2-4 h) - detect preformed enzymes
- Gas liquid chromatography (GLC)
  - Cell wall fatty acids
- MALDI-TOF MS
- 16S rRNA gene sequencing



# Special potency disk patterns

S ≥10 mm R <10 mm	Vancomycin 5 µg	Kanamycin 1000 µg	Colistin 10 µg
<i>B. fragilis</i> group	R	R	R
<i>Prevotella</i> spp.	R	R	V
<i>Porphyromonas</i> spp.	S	R	R
<i>Fusobacterium</i> spp.	R	S	S
<i>Veillonella</i> spp.	R	S	S


# Anaerobes of Medical Importance

# Pathogenic *Clostridium* and Assoc . Human Diseases

Species	Human Disease	Frequency
<i>C. difficile</i>	Abx assoc diarrhea, pseudom colitis	Common
<i>C. perfringens</i>	Soft tissue infections, food poisoning, enteritis necroticans, septicemia	Common
<i>C. septicum</i>	Gas gangrene, septicemia	Uncommon
<i>C. tertium</i> *	Opportunistic infections	Uncommon
<i>C. botulinum</i>	Botulism	Uncommon
<i>C. tetani</i>	Tetanus	Uncommon
<i>C. barati</i> , <i>C. butyricum</i>	Botulism	Rare
<i>C. histolyticum</i> *, <i>C. novyi</i> , <i>C. sordellii</i>	Gas gangrene	Rare
*aerotolerant		


# *Clostridioides (Clostridium) difficile*

- Asymptomatic colonization  
(~50% healthy neonates 1st yr)
- Disease is toxin-mediated \*\*All *C. diff* don't make toxin
  - Toxin A (enterotoxin)
  - Toxin B (cytotoxin)




## Disease:

- Antibiotic associated colitis
  - diarrhea
- Pseudomembranous colitis
- Toxic dilatation, perforation of colon and death



# *C. difficile* Toxin Detection

- Cell Culture Neutralization Assay
  - Add filtered stool + antitoxin to cells in culture
  - Toxin causes cytopathic rounding effect (CPE)
  - Antitoxin neutralizes CPE, showing specificity
- Expensive, 48 h TAT – not used routinely



# *Clostridioides difficile* Culture

- Toxigenic culture
  - Most sensitive (96%), long TAT
  - Yellow-white, ground glass colonies
  - Horse manure odor
  - Fluoresce chartreuse - UV light
  - Cytotoxin assay on 24 h broth cxs to show toxin is present



# Algorithms for *C. difficile* Detection

- Detection of *C. difficile* organism: GDH (glutamate dehydrogenase) Ag screen
  - If NEG, no further testing, saves \$\$ (>98% NPV)
  - If POS



- Toxin antigen assay
- Toxin assay, arbitrated by NAAT



- Antigen assays lower sensitivity if used w/o NAAT

# Algorithms for *C. difficile* Diagnosis

- Test for toxin by NAAT alone.
  - excellent sensitivity
  - poor specificity/PPV; over-diagnosis.
  - may detect colonization; non-viable bacilli
- NAAT, if positive do Toxin Ag assay
  - NAAT + Ag + : *C diff* present, making toxin
  - NAAT + Ag - : clinical correlation





# *C. difficile* Therapy

- Discontinue antibiotics, if possible
- Antibiotic therapy
  - Oral vancomycin
  - Fidaxomicin
- 10-20% relapse (persist spores, reinfect from environment)
- Fecal transplant
- Probiotics


# *Clostridium perfringens*

- Most common *Clostridium* spp. (20%) isolated in clinical specimens
- Common GI flora, transient skin flora
- Ubiquitous -intestinal tract animals, humans, type A can survive in soil
- Boxcar-shaped GPR
- Spores rarely seen



# *Clostridium perfringens*

- >16 toxins and enzymes
  - Types A-E based on production of  $\geq 1$  lethal toxin ( $\alpha, \beta, \epsilon, \iota$ )
  - Type A (alpha toxin) - most human infections
- Spectrum of disease
  - Soft tissue infections (cellulitis, suppurative myositis, myonecrosis or gas gangrene)
  - Food poisoning - heat-labile enterotoxin
  - Enteritis necroticans ( $\beta$  toxin)
  - Septicemia



# *C. perfringens* Food Poisoning

- Usually watery diarrhea with severe crampy abdominal pain.
- Uncommon: fever, chills, vomiting, headache.
- Incubation - 8-30 hours
- Associated with beef or chicken that has been stewed, roasted or boiled and allowed to stand and cool slowly.
- Type A enterotoxin – most common
- Self-limiting. Not treated.
- Stool not cultured in lab.



# Gas Gangrene (Myonecrosis)

- *Clostridium* spp.
  - *C. perfringens*, *C. histolyticum*, *C. septicum*, also *C. novyii*, *C. sordellii*, others
- Organisms contaminate wounds through trauma or surgery
  - Painful, foul-smelling, gas-filled, serosanguinous discharge.
  - Rapid necrosis
  - Shock, organ failure, 50% CFR if bacteremic
- Non-traumatic abdominal myonecrosis often due to *C. septicum*
  - Carcinoma, diverticulitis, surgery, chemo and other risk factors



Connor, Path Infect Dis, p. 524




# Alpha toxin of *Clostridium perfringens*

- Lecithinase activity
- Produced by all *C. perfringens*
- Lyses RBCs, plt, WBCs, endothelial cells
- Platelets & PMNs occlude blood vessels
- Increases vascular permeability




# *Clostridium perfringens*

- Rapid growth - detected after  $\leq 1$  day; divide every 8-10 min
- Nonmotile
- Double-zone  $\beta$ -hemolysis
  - Inner zone (complete): theta toxin
  - Outer zone (incomplete): alpha toxin (lecithinase)



# *Clostridium septicum*

- Gas gangrene, septicemia
- 70-85% w/ pos blood cx have underlying malignancy
- $\geq 68\%$  mortality
- Medusa-head colony morph on BBA 4-8 h; swarms by 24-48 h
- Ovoid (citron) subterm spores






# *Clostridium tetani*

- Tetanus
  - Puncture wounds
  - Tetanospasmin blocks release of neurotransmitters - inhibitory impulses are blocked in spinal cord & brain stem
  - Prolonged muscle spasms
  - Spasm of facial muscles
  - Spinal rigidity & spasm
  - Dx - clinical presentation
  - Rx - debride wound, MTZ, human tetanus IgG, vaccination



# *Clostridium tetani*

- Microscopy - poor sensitivity
- Prominent terminal spores (drumstick)
- Strict anaerobe - cx of wnd 30% sensitive
- Colony (slowly) spreading
- Tetanus toxin can be detected in serum



# *Clostridium botulinum*

- Botulism: life-threatening neuromuscular disease
  - also *C. baratii*, *C. butyricum*
  - 7 neurotoxins; A, B & E most common.
  - Prevents release of acetylcholine at neuromuscular junction
  - Blurred vision, impaired speech, flaccid paralysis
  - Foodborne, infant, wound
  - Home-canning; honey (infant)
  - **Public health reporting; bioterrorism concern.**
  - Dx: isolate org or detect toxin in food or stool, serum
  - Rx: vent, antitoxin, elim org (gastric lav, MTZ or Pen)

# *Actinomyces* spp.




- Gram pos branching rods; variable
- Colonies vary: white, grayish, red, pink, tan, yellow
- All facultative anaerobic except *A. meyeri*
- Actinomycosis, PID assoc w/ IUD, and pyogenic liver abscesses





# Actinomycosis

- Chronic granulomatous lesion forms abscess and draining sinus tracts
- Macroscopic sulfur granules
- Beaded, branching - usually in clinical specimens





# *Actinomyces* spp.

- *A. israelii* most common cause of actinomycoses in humans
- Smooth or molar-tooth colony
- Can be difficult to isolate - especially if specimen rich in microflora
- Also difficult to ID



# *Cutibacterium (Propionibacterium) spp.*

- Small GPR in short chains or clumps; diphtheroidal
- Catalase positive & indole pos = *C. acnes*
- *C. acnes*
  - Acne
  - Opportunistic infections assoc w/ prosthetic devices, IV lines






# *Bifidobacterium* and *Eubacterium* spp.

- Oropharynx, large intestine, vagina
- Low virulence
- Usually contaminants
- *Bifidobacterium*:  
pleomorphic GPR with  
bifurcated ends;
  - cocco-bacilli to long, branching rods
  - Colony is white to cream, convex, smooth, glistening



*Eubacterium* Gram Stain from thio broth;




*Bifidobacterium* (www.uaz.edu)




# *Lactobacillus* spp.

- Most facultative anaerobes; 20% obligate
- Some grow well aerobically
- Flora mouth, GI, GU tracts, foods
- Usually not clinically significant
- Rare cases of bacteremia, endocarditis




# *Lactobacillus* spp.

- Long and slender rods
- Often in chains
- Most nonmotile
- Catalase negative
- Colony alpha-hemolytic to gray, pinpoint to medium-sized.
- Many species are vancomycin resistant



# Gram-pos cocci: *Peptostreptococcus* et al

- Colonize oral cavity, GI tract, GU tract, skin
- Infection when spread to sterile sites
- Recently named genera
  - *Finegoldia magna*
    - Most pathogenic
  - *Parvimonas micra*
  - *Anaerococcus* spp.
  - *Peptoniphilus* spp.





- See fig 22.28 in text for simple ID scheme
- Vanc S, GLC, sequencing, MALDI-TOF MS for accurate ID

# Common Gram-negative Anaerobic Rods

- *Bacteroides fragilis* group (24 spp)
  - Most common
  - Bile resistant
  - More virulent and resistant than most other anaerobes
  - *Bacteroides* once included >50 spp
- *Parabacteroides*, other *Bacteroides* spp.
- Bile-sensitive species:
  - *Porphyromonas*
    - Pigmented, asaccharolytic; most indole pos
    - may fluoresce brick red
  - *Prevotella*
    - saccharolytic, bile sensitive
    - pigmented and nonpigmented
    - May fluoresce brick red

# *Bacteroides fragilis* group



- Grows well on BBE
- R to Van, Kan, Coli
- Circular, convex, grey-white colonies on BA
- Hydrolyze esculin (except *B. vulgatus*)





# *Bacteroides fragilis* group

- Pleomorphic, faintly staining gram-neg bacilli (may be vacuolized)





# *Campylobacter (Bacteroides) urealyticus*

- Pitting colony on BA
- Catalase neg
- Urease pos
- Bile sensitive; BBE neg
- Vanc R
- Kanamycin S
- Colistin S



# Pigmented *Prevotella* spp.- *Porphyromonas* spp. group

- GN coccobacilli
- (Brick-red) fluorescence under long-wave UV light
- Black pigmented colonies
- *Porphyromonas* spp.: vanc S  
no growth on LKV; spot  
indole +







# *Fusobacterium* spp.

- Disk
  - Vanc R
  - Kanamycin S
  - Colistin S
- *F. nucleatum* & *F. necrophorum* indole pos, nitrate neg
- *F. nucleatum*: thin GNR with pointed ends
- *F. necrophorum* & *F. mortiferum*: pleomorphic and swollen forms



# *Fusobacterium* spp.

- *F. nucleatum*
  - Common cause of abscess & BSI
- *F. necrophorum*
  - Lemierre's Disease – post-Strep infection. Septic thrombophlebitis, fever, chills, pain, swelling of throat
  - lipase pos on EYA




# *Fusobacterium* spp.

	Indole	Lipase	ONPG	BBE	Gram
<i>F. necrophorum</i>	+	+	-	-	Pleomorph
<i>F. nucleatum</i>	+	-	-	-	Slender, Pointed ends
<i>F. varium</i>	+	-	-	+	Regular shaped
<i>F. mortiferum</i>	-	-	+	+	Pleomorph swollen

# *Veillonella* spp.

- Small (<1.5  $\mu\text{m}$ ) GN diplococci
- Predominant anaerobe in human pharynx
- <1% of all anaerobes isolated
- Fluoresce red under UV, slow
- Disk
  - Vanc R
  - Kanamycin S
  - Colistin S



# Indications for susceptibility testing of anaerobes

- Serious infections: Brain abscess, endocarditis, osteomyelitis, joint infection, prosthetic device, bacteremia
- Normally sterile site
- Organisms considered highly virulent or unpredictable susceptibility profile
  - e.g., *B. fragilis* gp, *Prevotella*, *Fusobacterium*, *Clostridium*, *Bilophila*, *Sutterella* spp.
- Treatment failure (any anaerobe)

# M100, Appx D: Anaerobic Antibigram: 2013-16

*Bacteroides fragilis* group (7 spp.), n=2580.

% Susceptible, depends on species

- Ampicillin-sulbactam, 45-84%
- Piperacillin-tazobactam, 87-96%
- Cefoxitin, 13-100%
- Ertapenem, 82-84%
- Imipenem, 97-100%
- Meropenem, 93-100%
- Metronidazole, 100%
- Clindamycin, 26-53%
- Moxifloxacin 31-62%

# CLSI M100, Appx D: Anaerobic Antibigram: 2013-16

- Other anaerobes; %Susceptible

[*Prevotella*, ANA GPC, *Fusobacterium*, *Cutibacterium*, *Clostridium* spp., not *C. diff*]

- Ampicillin-sulbactam, 97-100%
- Piperacillin-tazobactam, 94-100%
- Imipenem, 94-100%
- Meropenem, 98-100%
- Metronidazole, 95-100% (*C. acnes* 0% Susc)
- Clindamycin, 53-97%
- Moxifloxacin, 62-95%


Anaerobic gram-positive cocci include *Peptococcus*, *Peptostreptococcus*, *Finegoldia*, *Peptoniphilus*, and *Anaerococcus* sp

# *C. acnes* should be susceptible to Penicillin

Benzylpenicillin / *Propionibacterium acnes*  
International MIC Distribution - Reference Database 2017-10-01

**EUCAST penicillin clinical breakpoints for anaerobes: S  $\leq 0.25$ , R  $> 0.5$  mg/L**

**CLSI penicillin breakpoints for anaerobes: S  $\leq 0.5$ , I = 1, R  $\geq 2$   $\mu\text{g/mL}$**




European Committee on Antimicrobial Susceptibility Testing. Data from the EUCAST MIC distribution website, last accessed 10/1/2017 <http://www.eucast.org>



# Acknowledgement

- Dr. Sandra Richter

# Toxicity of Oxygen



# *Clostridium difficile*

- Most common cause of healthcare associated diarrhea
  - Doubled incidence since 1996; outbreaks common
  - Cost >\$1 billion/yr
  - 027/NAP1 strain quinolone R outbreak strain
  - Risk factors: Clinda, amp, amox, ceph, quinolones (up to 8 wks later)
- Most toxigenic isolates A+/B+
  - Enterotoxin (toxin A)
  - Cytotoxin (toxin B)