


LAB BASICS

RYAN COLLISON, MLS (ASCP)^{MLS}


TOPICS

- [Specimen Processing](#)
- [Laboratory Supplies](#)
- [Pipettes](#)
- [Reagents](#)

















SPECIMEN IDENTIFICATION

- Throughout process specimen labeling is vital
 - 2 unique IDs
 - First and last name, birth date, MRN, Contact
- Making aliquots
 - Secondary container is clean
 - Pipettes clean
 - Good technique
 - Transfer one specimen at a time
 - Visual signals (physically moving tube etc)



SPECIMEN CONTAINERS

- Phlebotomy tubes
 - Whole blood
 - Anticoagulated, kept homogenous or re-mixed
 - Plasma
 - Anticoagulated and spun in centrifuge
 - Citrate, heparin, EDTA, K-oxalate
 - Serum
 - Whole blood clots, spun in centrifuge
- Functional additives
 - NaF, preservatives

Citrate		Sodium citrate, plastic	
		Sodium citrate & CTAD, glass	
		Sodium citrate, ESR, glass	
		ACD	
Serum		Serum with clot activator (silica particles)	
		Serum thrombin	
		BD RST (serum with gel)	
		BD SST™ II Advance (serum with gel)	
Heparin		Lithium & sodium heparin	
		BD PST™ II (plasma with gel)	
Haematology		EDTA	
Crossmatch		Blood banking	
PPT		BD PPT™ EDTA with gel	
Trace elements		Trace elements with EDTA	
		Trace elements serum with clot activator (silica particles)	
Glucose		Glucose	

SPECIMEN PROCESSING

- Centrifugation
 - Separate particles from solution
 - Concentrate cellular components
 - Enhance Ag-Ab reactions
 - Precipitate protein components
 - Separate liquid phases based on density



CENTRIFUGATION

- Rotor
 - Holds specimen
- Motor
 - Speed regulated by voltage input
 - Tachometer measures, adjusts speed
- Drive shaft
 - Connects motor and rotor
- Timer, Brake
 - Reaches end, slows down rotor



CENTRIFUGATION

- Centrifugal Forces
 - Speed of centrifuge (RPM) does not tell us force
 - Relative Centrifugal Force (RCF) $RCF = 1.12 r \times \left(\frac{RPM}{1000}\right)^2$
 - Radius in mm
 - Smaller radius at same sp
 - Resulting unit is x g (multiples of g-force)
 - Resulting force is high, centrifuge **MUST** be well-balanced



LABORATORY SUPPLIES

- Glassware
 - Borosilicate
 - Aluminosilicate
- Beaker
- Flask
 - Erlenmeyer
 - Volumetric
- Graduated Cylinder




- Pipettes
- Burets
- Plastic
 - The new glass!



LABORATORY SUPPLIES

HOW TO USE: BEAKERS

- Holding bulk volume of solutions
 - If using small volumes consider using a medicine cup
- Waste containers
- All glassware can be labeled with sharpie
 - Alcohol prep pads clean off sharpie
- Not used for measurement



LABORATORY SUPPLIES

HOW TO USE: FLASKS

- Erlenmeyer
 - Use similar to beakers
 - If sealed by parafilm may mix by swirling
 - Not used for measurement



- Volumetric
 - Holds a single exact volume
 - One calibration mark and tolerance listed
 - Application? i.e. make a 10 g/L NaCl
 - In 1L flask put 10 g NaCl
 - Add almost 1 L of water, say 950 mL
 - Add rest of water dropwise with transfer pipette to meniscus
 - Insert stopper, invert and swirl to mix



LABORATORY SUPPLIES

HOW TO USE: GRADUATED CYLINDERS

- Use class A for measurements
 - These are very inexact
 - Convenient for pouring
 - Filling 50 mL volumetric flask? Pour 45 mL using graduated cylinder
- Short answer?
 - Don't use for important things



PIPETTES

- Class A
 - NIST standards
 - Calibrated
 - CAP requires using calibrated pipettes
- TC versus TD
 - TD: What you think of when you think pipette
 - TC: Requires washing into solution



PIPETTES

- Transfer- transfers single known volume
 - Volumetric pipettes
 - Ostwald-Folin ||
- Measuring- labeled in increments
 - Serological ||
 - Mohr
- Micropipette
 - Adjustable or fixed volumes ≤ 1 mL


Volumetric



Serological



Mohr



Ostwald-Folin




VOLUMETRIC PIPETTE

- Elongated bulb
- Tapered Tip
- Single Calibration mark
 - Dispenses ONLY the volume stated
- Self draining
- Most precise & accurate pipette
 - Suitable for calibrators



SEROLOGICAL PIPETTE

- Uniform width
- Wider tip
- Measuring pipette
 - Calibrated through the tip
 - Can dispense any volume for which a graduation exists
- Usually has etched ring (or 2) at top of pipette
 - blowout



MOHR PIPETTE

- Similar appearance to serological
- NOTE GRADATIONS
 - Stopping before tip
- Never blowout, never drain to tip
- Why Mohr?



OSTWALD-FOLIN PIPETTE

- Appearance of volumetric pipette
 - Bulb closer to tip
- Frosted rings
 - Blowout
- Why?



PIPETTING PROCEDURE

1. Inspect pipette
 1. Chips, cracks, dirt
2. Using bulb draw liquid up past calibration mark
3. Wipe pipette with kimwipe
4. Holding vertically: allow liquid to drain down to calibration mark (QS'ing)
 1. Place finger to hold meniscus on calibration mark
 2. Move to receiving vessel
5. Allow to drain unrestricted to bottom – IF draining to another calibration mark, stop here
6. Touch tip to side OR blowout depending on pipette

PIPETTING RULES

- Do Not Mouth Pipette
 - Why do we still have to say this? I don't know
- Always check the pipette for chips/cracks/dirt
- Lay pipettes on clean surface
- Make sure there are no bubbles or beads of liquid when using
- Place in pipette soak tip down
- NEVER out of a primary container
 - Use a labeled pourover



AUTOMATIC PIPETTE

- 2 flavors- adjustable and fixed volume
 - Adjustable should be used at top of adjustable range if possible
- Holding pipette vertically depress plunger to first stop
- Place tip in liquid
- Slowly release plunger to top
- Raise pipette, inspect tip for bubbles, wipe outside of tip (DO NOT BLOT ORIFICE)
- In receiving container touch tip to side, depress plunger steadily
 - If second stop, hesitate then depress to second stop.
- Remove pipette, eject tip, release plunger



REAGENTS

- Purity of Inorganics
 - Ultrapure- chromatography, atomic absorption, molecular diagnostics
 - Analytical reagent- Amer. Chem. Soc. Definition actual impurities are listed with %
 - USP/NF- For pharmaceutical use, it won't kill anyone but we don't know purity
 - Chemically pure, pure grade- starting point for manufacture, do not use
- Purity of organics
 - Chromatographic grade- Depends on refinement level
 - Reagent grade- certified to have < % of contaminants
 - Spectroscopic- spectroscopically pure
 - CP- almost reagent grade



HANDLING REAGENTS

- Store properly
 - RT = 21°C
- Never sample directly from bottle
 - Pour into cup, beaker, or boat
 - Dispose of according to SDS
 - Expiration dates MUST be honored

**EXPIRATION
DATE: --/--**

WATER PURITY



- Water Water Everywhere, but a not a drop for quantitative clinical chemistry (amiright?)
 - Ultrapure- used for HPLC, pharmaceuticals
 - Type I- Chemistry procedures, cell culture, electrophoresis, reagent/control prep, final rinse of glassware
 - Type II- qualitative chemistry, most heme micro and immunology procedures
 - Type III- used to make Type I, II water, used in first rinses of glassware
- Purity mainly assessed using conductivity

