


Electrophoresis testing

Sonja Bruketa, MLS(ASCP)

Principles of electrophoresis

- Electrophoresis testing uses an electrical field to separate molecules based on their charge





Oncohemat Key

Electrophoresis Techniques | Oncohemat Key

Factors which effect electrophoretic mobility

- 3 major components of electrophoresis techniques
 - Support media (paper, cellulose acetate, agarose) – usually net negative charge that doesn't move
 - Buffer/solvent - usually net positive charge that flows toward cathode (electro-endosmosis)
 - Patient sample – most major serum proteins net negative charge at pH used in assays (8.6-10.0)




- Migration on gel depends on strength of negative charge for a given protein vs the endosmotic “current” they’re swimming against

Clinical applications of electrophoresis


- 3 main clinical scenarios where ELP testing is used
 - Identification of monoclonal immunoglobulins (multiple myeloma, MGUS, amyloid)
 - Identification of oligoclonal bands in CSF (multiple sclerosis, CNS inflammatory disease)
 - Identifying the source of an abnormal serum enzyme level (ALP, CK, LDH, A1AT)
- 2 major techniques used for this testing
 - Protein electrophoresis: separate by charge, visualize with non-specific stain
 - Immunofixation: separate by charge, visualize with antibody-specific fixation and staining
- Testing is usually done on serum or urine, although other types may be used in specific situations

Protein electrophoresis

- Serum/urine applied to agarose gel and current applied.
- Electrophoresed sample stained with non-specific (acid blue) protein stain
- Banding pattern is examined to look for abnormalities relative to normal distribution





Normal serum protein migration pattern




- The locations and intensity of bands are fairly consistent between individuals, allowing identification of abnormal bands

ELP patterns in multiple myeloma



Advantages	Limitations
High-throughput, inexpensive screen	Abnormal PE bands can be due to other causes than monoclonal gammopathy
Provides semi-quant estimate of amount of M protein present	M proteins can "hide" behind normal serum proteins
Can detect other abnormalities besides monoclonal gammopathy	Technique is less sensitive than IFE for low-level M proteins

ELP patterns in other common diseases



Low albumin/ $\alpha 2$ due to decreased liver synthetic function

Beta-gamma bridging due to increased IgA

Low albumin due to renal losses

$\alpha 2$ increase from liver compensatory increases in synthesis

$\alpha 2$ increase due to hepatic synthesis of acute phase proteins


$\alpha 2$ increase due to hepatic synthesis of acute phase proteins

Polyclonal increase of immunoglobulins

- However, without an IFE one cannot definitively rule out a monoclonal Ig as the cause of a band


Common serum EP artifacts

Plasma




Sharp band in γ from fibrinogen

Hemolysis



Shoulder on α_2 region from hemoglobin-haptoglobin complex

Hepatitis





Multiple faint γ bands on polyclonal background from chronic inflammation

V8 Nexus Capillary Electrophoresis


**Fully automated, faster,
more connected and
flexible**

**Patient's results can be
integrated with gel results
in tracking of treatment**





Immunofixation electrophoresis

- Serum/urine applied to agarose gel and current applied.
- Antisera applied to gel which precipitates the appropriate isotype (G, A, M, K, or L)
- Unreacted proteins washed away and precipitated complexes stained with non-specific (acid violet) stain




Normal and abnormal IFE patterns




Advantages	Limitations
Confirms that abnormal band is a monoclonal immunoglobulin	Does not quantitate amount of M protein present
Detection not affected by normal serum protein bands	Can detect low-level bands on unclear significance
Provides information on isotype	Lower throughput

Isoelectric focusing

- Paired serum and CSF are applied to agarose gel and current applied.
- Samples migrate through a pH gradient where protein charges change
- Proteins stop migrating when they reach zero charge (isoelectric point)
- Provides improved separation of proteins with similar charges at biological pH




IEF in multiple sclerosis





- Normal appearance is no bands in CSF or serum
- Bands in CSF > serum :
intrathecal synthesis of IgG
 - Can be due to causes other than MS
- Bands in CSF = serum : usually systemic inflammatory process with crossover into CSF
 - Crossover can be *in vivo* or due to contamination during spinal tap

Isoenzyme testing (alkaline phosphatase)



Alkaline phosphatase isoenzyme testing

- Serum is treated with neuraminidase (cleaves sialic acid to improve bone/liver separation)
- Samples electrophoresed on agarose gel in constant pH buffer
- Gels treated with NBT (nitro blue tetrazolium) which is cleaved to blue stain by alkaline phosphatase
- Fractions in bone, liver, and other are quantitated as % of whole



Conclusions

- Electrophoresis techniques take advantage of charge differences to characterize proteins
- Immunofixation can be used to increase the sensitivity and specificity of testing for monoclonal immunoglobulins
- Although most testing can be done with routine techniques, specialized adaptations such as isoelectric focusing and isoenzyme analysis can be used for specific clinical conditions.